

Point of Care Testing Quiz

1. Lancets can be used on multiple residents/patients.
 - a. True
 - b. **False**

2. Only shared glucometers have to be cleaned and disinfected after each.
 - a. True
 - b. **False**

3. Glucometers should be cleaned and disinfected after each use with:
 - a. An alcohol wipe.
 - b. **A disinfecting wipe as specified by the IFUs.**
 - c. Glucometers do not need to be cleaned and disinfected after each use.
 - d. A disinfecting product that your facility has available.

4. Used lancets should be discarded:
 - a. In the closest garbage bin
 - b. In a glove that was used during the procedure
 - c. **In a sharps bin**
 - d. None of the above

5. All unused supplies that were not needed for the procedure but were taken into the resident/patient room should be:
 - a. **Discarded.**
 - b. Taken back to the clean supply storage area.
 - c. Taken back to the medication cart and stored.
 - d. Used on the next resident/patient.